





# Tips & Tricks



# What is today's objective?



- to introduce possibilities of Mnova@IOCB
- to simplify your Mnova workflow
- to show some features of Mnova program you might not be aware of
- to find your demand for Mnova workshops



866 pages

MestReNova Manual

© 2021 MESTRELAB RESEARCH Last Revision: 3-May-2021

#### https://www.youtube.com/Mestrelab



tutorials, webinars, tips





# Where to download?

https://mestrelab.com/download/mnova/

Licenses Samba server: 147.231.122.7/nmruser/MestrecLic

- IOCB holds unlimited number of licenses for NMR, MS and UV/Vis plugins
- You can install Mnova on your laptop works 90 days outside IOCB

# Available plugins/licenses

# $\rightarrow$ Help/License Manager

Li	icense I	Manager							?
Hos	st ID:								
		L5NV1-1MHH123W-V	D20X-MA0MLZBK						
lice	nses								
	State	Plug-in	Issued By	Licensed To	Туре	Issued Date	Days to Expire	Update Days	Valid Day
1	8	EIViS	Mestrelab Research S.L.	UOCHB AV CR	campus	pá kvě 14 2021	Never	973	N/A
2	0	Mass	Mestrelab Research S.L.	UOCHB AV CR	campus	pá kvě 14 2021	Never	973	61
3	0	Mestrelab Predictor	Mestrelab Research S.L.	UOCHB AV CR	campus	pá kvě 14 2021	Never	973	N/A
4	0	Mnova qNMR	Mestrelab Research S.L.	UOCHB AV CR	campus	pá kvě 14 2021	Never	973	61
5	0	NMR	Mestrelab Research S.L.	UOCHB AV CR	campus	pá kvě 14 2021	Never	973	61
6	0	NMRPredict Desktop	Mestrelab Research S.L.	UOCHB AV CR	campus	pá kvě 14 2021	Never	973	61
7	0	Random Forest Predictor	Mestrelab Research S.L.	UOCHB AV CR	campus	pá kvě 14 2021	Never	973	61
8	0	Reaction Monitoring	Mestrelab Research S.L.	UOCHB AV CR	campus	pá kvě 14 2021	Never	973	61



# Edit/Preferences

General: large icons, vertical tabs NMR: import parameters (e.g. appodization, zero filling) *View/Toolbars View/Tables* – parameters, compounds, multiplets, peaks *View/Panels* – data browser, cursor info *View/Pages* Ctrl+F2 Shortcuts



# to see them go to: *Help/Shorcuts*

#### Important shortcuts for NMR:

- I manual integration
- J manual multiplet analysis
- K manual threshold for peak pickingCtrl+K peak picking peak by peak
- L reference
- **R** graphic reference
- **W** apodization
- **B** baseline correction
- Shift+P phase correction
- **C** crosshair
- **E** expansion
- decrease intensity
- + increase intensity
- F full spectrum
- Z zoom

Sho	tcuts		
	Command	Shortcut	
1	🧩 Analysis > Assignments > Manual Assignment	Α	
2	✿ Analysis > Assignments > Swap Assignments	s	
3	Analysis > Integration > Delete Manually	Ctrl+Shift+I	
4	Manual Manual Manual	I	
5	₩ Analysis > Multiplet Analysis > Manual	J	
6	Analysis > Peak Picking > Delete Manually	Ctrl+Shift+K	
7	🞽 Analysis > Peak Picking > Manual Threshold	к	
8	📌 Analysis > Peak Picking > Peak by Peak	Ctrl+K	
9	★ Analysis > Reference > Graphic Reference	R	
10	★ Analysis > Reference > Reference	L	
11	Documents > Close	Ctrl+F4	
12	T Edit > Annotate > Text	т	
13	Edit > Copy	Ctrl+C	
14	📑 Edit > Create New Page	Ctrl+M	
15	🔏 Edit > Cut	Ctrl+X	
16	🗱 Edit > Delete	Del	
17	🖺 Edit > Paste	Ctrl+V	
18	O Edit > Preferences	Ctrl+,	
19	C <sup>★</sup> Edit > Redo	Ctrl+Y	
20	Edit > Select All	Ctrl+A	
21	${\underset{A_{A}}{\overset{\bullet}}} \  \  \  \  \  \  \  \  \  \  \  \  $	Ctrl+Shift+-	
22	${}_{\scriptscriptstyle A}{}^{{}_{\scriptscriptstyle A}}$ Edit > Text Format > Increase Font Size	Ctrl+Shift++	
23	★ Edit > Undo	Ctrl+Z	
24	🔀 File > Export to PDF	Ctrl+D	
25	File > New	Ctrl+N	
26	階 File > Open Directory	Ctrl+I	

	Command	Shortcut	
	Command	Shortcut	
27	File > Open	Ctrl+O	
28	File > Print	Ctrl+P	
29	File > Save	Ctrl+S	
30	Help > Contents	F1	
31	Integral Manager	Shift+I	
32	Multiplet Manager	Shift+J	
33	≽ Processing > Apodization	w	
34	▶ Processing > Baseline > Baseline Correction	В	
35	脊 Processing > Fourier Transform	Shift+F	
36	↑ Processing > Phase Correction > Manual Correction	Shift+P	
37		с	
38	<u>ී</u> ෆී View > Cuts > Cut	х	
39		v	
40	📇 View > Expansion	E	
41	View > Full Screen	F11	
42	✓ View > Intensity > Decrease	-	
43	M View > Intensity > Fit to Highest Intensity	н	
44		+	
45	View > Pages	Ctrl+F2	
46	🗛 View > Pan	Р	
47	View > Zoom > Full Spectrum	F	
48	💘 View > Zoom > Manual Zoom	М	
49	📌 View > Zoom > Next Zoom	Shift+Right	
50	View > Zoom > Previous Zoom	Shift+Left	
51	😻 View > Zoom > Zoom In	Z	
52	View > Zoom > Zoom Out	Shift+Z	



#### **Dataset opening:**

- 1. File/Open (Ctrl+O) or File/OpenDirectory (Ctrl+I)
- 2. Drag and drop the folder with NMR data (the folder can contain structure in \*.mol file will be open automatically)
- 3. Data browser panel selection of location containing data

**Customizing the spectrum appearance** (graphical properties) – *Edit/Properties* or right mouse button click on spectrum and then *Properties* 

### Essential steps in 1D NMR data processing (in *Processing*):

- 1. Fourier transformation
- 2. Zero filling improves digital resolution
- 3. Appodization (W) enhances sensitivity (exponential function) or resolution (exponential function + Gaussian function); in order to see the FID+function click on *View/FullView*
- 4. Phase correction (Shift+P)
- 5. Baseline correction (B), multipoint baseline correction
- 6. Referencing (L)
- 7. Peak picking (K setting threshold, Ctrl+K peak by peak), *PeakPicking/Options*
- 8. Integration (I)



- 1. use Data Browser for storing of your NMR data
- 2. create your favorite graphical properties of spectra and save it
- 3. create your favorite **Processing Template** for each type of spectrum (<sup>1</sup>H, <sup>13</sup>C, <sup>31</sup>P...) and save it

# Other "tricks":

- several types of spectra on single page
- expanding of spectrum View/Expansion or E
- finding experiment parameters *View/Tables/Parameters*
- chemical structure (supported format of chemical structures: \*.mol \*.cdx \*.sdf)
- signal suppression *Processing/Signal Suppression*



# MestReNova enables easy creation of Supporting Info:

- 1. Pasting a spectrum into MSWord- possibility of later formatting the spectrum
- 2. Export to PDF *File/Export to PDF* or **Ctrl+D**
- 3. Copy as a picture Edit/Copy Special/Copy as Image

# Peaks report generation:

- 1. Creating of peak list by peak picking
- 2. Peak report generation *Scripts/Report/Peaks* or in Peaks Table *Report Peaks*
- 3. Changing style of peak report in Peaks Table Setup Report

# Multiplets report generation:

- 1. Creating of multiplet list by multiplet analysis (J)
- 2. Peak report generation *Scripts/Report/Multiplets* or in Multiplets Table *Report Multiplets*
- 3. Changing style of multiplet report Multiplets Table Setup Report



- 1. select spectra in Pages Table
- 2. Go to Stack/Stack Spectra
- 3. To work with stacked spectra  $\rightarrow$  activate toolbar *Stacked*



Graphical properties of the stacked spectra  $\rightarrow$  *Properties/Stacked* 

# Reaction kinetics monitoring



# In stacked mode: Advanced/Data Analysis/Create/Integrals Graph



- Paste molecular structure or drag and drop file with your structure (\*.mol \*.cdx \*.sdf) to *Compounds* Table and run prediction of <sup>1</sup>H, <sup>13</sup>C, and other nuclei *Predict/1H Spectrum* or <sup>13</sup>C Spectrum
- You can try also for other nuclei (<sup>11</sup>B, <sup>15</sup>N, <sup>19</sup>F, <sup>29</sup>Si, and <sup>31</sup>P)
- Predict/Prediction Options selection of e.g. solvent, frequency, number of datapoints....
- Compare predicted spectrum with your experimental spectrum



• paste published spectrum into *Multiplet Report to Spectrum* script and generate spectrum (*Scripts/NMR Tools/Multiplet Report to Spectrum*)

example for JM040 – 2,6-diamino-2'-deoxyadenosine

**2,6-Diamino-2'-deoxyadenosine (2)**—Supplied by Raylo Chem. Edmonton, Alberta, Canada. <sup>1</sup>H NMR (300 MHz,  $d_6$ -DMSO):  $\delta$  2.17 (m, J = 2.7, 6.0, -13.2 Hz, 1H), 2.60 (m, J = 5.7, 8.4, -13.5 Hz, 1H), 3.52 (m, 1H), 3.59 (m, 1H), 3.84 (m, 1H), 4.35 (m, 1H), 5.26 (m, 5'-OH, 1H), 5.29 (m, 3'-OH, 1H), 5.76 (br s, 2 NH<sub>2</sub>, 2H), 6.17 (dd, J = 6.0, 8.1 Hz, 1H), 6.76 (br s, 6 NH<sub>2</sub>, 2H), 7.92 (s, 1H), <sup>13</sup>C NMR (75 MHz,  $d_6$ -DMSO):  $\delta$  39.31, 61.91, 70.94, 83.03, 87.57, 113.39, 135.69, 151.12, 156.08, 159.94. MS: m/e 267 (M + 1)<sup>+</sup>. UV max: 255, 280 nm.

Journal of Pharmaceutical Sciences 1994, 83 (4), 525-531.

• compare with experimental spectrum



**NMR** SPECTROSCOPY offers Mnova workshops covering various topics:

- 1D NMR processing
- 2D NMR processing
- multiplet analysis
- spectra assignment
- presenting NMR data
- reaction monitoring and qNMR
- spin system simulation
- Mnova scripts
- troubleshooting

#### workshops concept:

- 10 attendees with laptops
- topics will be specified upon demand
- hands-on solving of particular problem

