

August 2017

Standard Operating Procedure

Optical setup of Aria IIu

Pages 1/1

	<i>Mirror</i>	<i>Filter</i>	<i>Parameter</i>	<i>Fluorochromes</i>
Violet 405nm	502LP	530/30	Alexa Fluor 430	AF430
	-	450/40	DAPI	DAPI, Pacific Blue, AF405
Blue 488nm	7735LP	780/60	PE-Cy7	PE-Cy7
	635LP	695/40	PerCP-Cy5.5	PerCP-Cy5.5
	595LP	610/20	PE-texas red	PI, RFP, PE-Texas red
	556LP	575/26	PE	PE
	502LP	530/30	FITC	FITC
	-	488/10	SSC	
Red 640nm	735LP	780/60	APC-Cy7	APC-Cy7, APC-H7, AlexaFluor 750, Zombie NIR
	-	660/20	APC	APC