

Cell Request Form

IOCB PRAGUE

Group of Biochemical Pharmacology

Contacts: Helena Mertlíková Kaiserová (organization issues, primary cultures, other specific requests)
Karolína Müllerová (requests on the layout of requested cells)
Alžbeta Magdolenová (PBMC)

Location of Tissue culture laboratory: Room A 2.28

Tel: Helena Mertlíková Kaiserová 220-183-114, office A 2.31
Karolína Müllerová 220-183-521, office A 2.99
Alžbeta Magdolenová 220-183-528, office A 2.30

E-Mail: **cellservice@uochb.cas.cz**

Contact detail:

Your name, tel. no and e-mail address:
Name of the lab and group leader:
Short description of the project (e.g. cytotoxicity, microscopy...) :

Cells to be provided:

I do not have specific cell line in my mind, please help me to choose based on the following criteria:

Species:	Human	Mouse	Rat	
Adhesiveness:	Adherent	Suspension		
Tissue:	Blood cells	Lung	Heart	Liver
	Breast	Cervix	Brain	Kidney
	Intestine	Prostate	Connective tissue	

I wish to get specific cell line (check immediate availability beforehand!):

CCRF-CEM	HL-60	HepG2	HeLa	T47D
U2OS	U87-MG	MCF	DU-145	LoVo
Caco-2	A549	J774A.1	RAW264	THP-1
H9c2	MOLT-4	1301	Jurkat	HEK-293
L1210	SH-SY5Y	CA-46		

Seeding density (number of cells per well/ per cm² / per ml of suspension):

Media volume per well (ul):

Special media requirements: Yes No

If yes, please specify:

Date when the cells should be ready*:

Requested Plasticware:

Multiwell plate					
Number of wells:	6	12	24	96	384
Plate Color:	Transparent		Black	White	
Optical bottom (plastic**):	Yes		No		
Petri dish	6-cm	10-cm	Other (please specify)		
Flask	T25	T75	T150		
Tube	1,5 ml microtube (Eppendorf)				
	15 ml Falcon		50 ml Falcon		

*pick-up in the afternoon hours unless otherwise agreed on

**if you have highest demands on optical quality of the plate (high resolution microscopy), glass-bottom plates are recommended. Please note that most cell types will not grow on glass surface and would require additional coating with matrices like poly-D-lysine, collagen, gelatin, laminin etc. If you are unsure with the choice of the right plate for your application, we are ready to help.